



# RNA Sequencing: Improving Diagnostic and Therapeutic Outcomes for Genetic Testing

W.A. Langley, Y. Cuperus, N.A. Rouse, Jiyong Wang, P.L. Nagy


## Introduction

- Next generation sequencing of DNA is extremely powerful in identifying the cause of genetic disorders. Finding a rare, "in silico" disruptive variant is not sufficient to establish its pathogenicity. Reporting of variants of unknown significance (VUS) causes anxiety and increases the workload of the care providers. Functional tests, such as transcriptome analysis/ RNA sequencing can help in the evaluation of the molecular consequences of both SNVs and CNV.
- RNA-sequencing<sup>1</sup> can demonstrate the pathogenic nature of deletions, frameshifts or non-canonical splice-site variants by revealing the effects of these changes on the level and content of transcripts present in the affected tissue. It is even useful in the diagnoses of imprinting disorders, irrespective of the underlying molecular mechanism. If the gene is widely expressed, blood or other easy to sample tissue types can be used with similar efficiency. Transcriptome analysis also helps to identify patients who could benefit from disease therapies based on splice-modulating therapeutic approaches.
- We have developed validated and have been using clinical transcriptome sequencing to assist with the diagnosis of constitutional disorders since last November. In our experience transcriptome analysis can rule in or rule out the proposed splicing disruptive effect of VUS's in almost all cases where appropriate tissue is available. In addition, in some cases it can also provide diagnosis in the absence or instead of DNA-based diagnostic approaches (Prader Willi syndrome and SMA testing).

## Visualization with Integrated Genome Viewer (IGV)



Fig 1: DMD diagnosis



Prader Willi Syndrome

- Prader Willi critical region on 15q11.2-q13
- Abnormal parent-specific imprinting, 3 known mechanisms: paternal deletion, uniparental disomy 15 (UPD), and imprinting defect (ID)
- Severe hypotonia and feeding difficulties early in life
- Later excessive eating and central obesity
- Delay in motor milestones and language development
- Cognitive impairment, with behavior problems
- Short stature
- Hypogonadism, incomplete pubertal dev, and infertility

### Prader-Willi syndrome : Genetic mechanisms





Fig 4: Assessment of splice sites VUS



## RNA Sequencing Test Options

### MNG Transcriptome™

Complete RNA transcriptome sequencing that can be ordered at the time of initial testing as a reflex to any NGS panel, or in conjunction with our MNG Exome™.

### Panel-Specific RNA Sequencing

RNA sequencing analysis can be ordered as an add-on to any of MNG Laboratories' NGS panels and will only include analysis of the genes listed in the panel.

### Gene-Specific RNA Sequencing

Analysis of up to five genes for variant investigation following a report with a VUS identified in a splice site or intronic region, or the variant is noted to have a potential effect on splicing. Previous sequencing analysis can be done at any laboratory.

## Why Perform Transcriptome Testing - Summary

- Transcriptome is the summary of all different genetic changes
  - Regulatory mutations
    - Enhancer and promoter mutations
  - Large deletions of coding regions
    - Transcript level is decreased by half
    - In imprinted loci transcript is absent
  - Intragenic indels
    - Alternative initiation and termination sites
    - Exon skipping
    - Nonsense mediated decay
  - Intragenic SNVs
    - Splice site efficiency modulation
    - Premature stops activating nonsense mediated decay
- With the development of splicing modifiers, it is important to assess the exact processing outcome of a specific variant

## Methodology - Design

- Next generation sequencing library preparation by TruSeq Stranded Total RNA Library Prep Kit/ Ribo-Zero Human (Illumina)
- Validated for blood, muscle, brain/nerve, fibroblast culture
- RNA sequencing analysis is run on a HiSeq or NovaSeq instruments (Illumina)
- Housekeeping genes were determined as the set of genes that had average expression and low variance across all tissue types
- Reportable genes were defined as genes that have expression > FPKM1 and low variance in a given tissue
- Reference ranges were established for junction, gene, and transcript expression using the tissue specific control sets
- To determine the preferred sample type, the human protein atlas<sup>3</sup> is used as one of the sources to look at gene expression in different tissues

## Determining Tissue Type for Testing

MNG uses multiple criteria and a team approach to help determine the best sample type to submit based on clinical relevance to the patient and variants

Internally developed RNA expression profiles and reportable ranges across all genes and tissue types

Publicly available databases for protein and RNA expression

Clinical information relevant to the patient's phenotype and symptomology

## Quality control of RNAseq samples

- High quality reads: > 85 % of all reads Q30 or greater
- Total yield per sample: > 10 Gb
- Alignment statistics
  - Alignment rate: > 90 %
  - Number of uniquely mapped reads: > 50 million
  - rRNA contamination percentage: < 2 % of all mapped reads
- Other metrics tracked
  - Transcript reconstruction statistics
  - Exon-exon junction reconstruction statistics

Fig 2: Prader-Willi diagnosis





Fig 3: SMA diagnosis



## Tissue Specific Shipping Conditions

Tissue Types	Blood	Fibroblasts and Lymphoid Cell Lines	Nerve, Brain, Muscle Tissue
<b>Sample Preparation and Collection</b>	MNG Provided Tempus™ Blood RNA Collection Tube (Thermo Scientific)	Cultured, T-25 flask Preferred: cell pellet spun down and frozen	50-75 milligrams of fresh, frozen tissue
<b>Shipping Conditions</b>	Room temperature, overnight shipping	5-7 pounds dry ice, overnight shipping	5-7 pounds dry ice, overnight shipping

Please call MNG Labs at 678-225-0222 for any questions or assistance

## Bioinformatics

- Tools used:
  - HISAT (hierarchical indexing for spliced alignment of transcripts), StringTie and Ballgown
  - <http://ccb.jhu.edu/software.shtml>
    - Aligns reads to a genome
    - Assembles transcripts including novel splice variants
    - Compute the abundance of these transcripts and splice junctions in each sample
    - Compare experiments to identify differentially expressed genes and transcripts

## References

- Cummings et al, Improving genetic diagnosis in Mendelian disease with transcriptome sequencing. Sci Transl Med, 2017; 9(386) PMID 28424332
- Pertea et al, Transcript-level expression analysis of RNA-seq experiments with HISAT, StringTie and Ballgown. Nat Protoc, 2016; 11(9): 1650-1667 PMID 27560171
- www.proteinatlas.org